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Statement:

The method of Tengzhou City Xiangying Analysis Technology Co., Ltd. (Chromatographic engineering research and analysis and test center) exclusive research and development with the nature of the patent, protected by the patent law, where the identical, belong to the fake, Xiangying technology will be prosecuted according to law the legal liability of the infringer.

preface

Years of grain and oil testing work enable us to clearly understood: strict implementation of national standards, in accordance with the provisions of the standard procedures and steps were detected, is the bounden duty of every inspection of grain and oil workers. However, in accordance with the people's Republic of China national standard GB5009.37 the edible method for analysis of hygienic standard of plant oil, residual solvent determination steps are performed to detect, often appear some non-human but inevitable problems. Therefore, Tengzhou City Xiangying Analysis Technology Co., Ltd. according to all the year round in the actual testing process encountered some problems, of GB were some amendment and supplement, developed a new approach to the analysis, and achieved good results, special to share with the new and old customers. Deficiencies in the national standard:

1. GB successful sample preparation using N, N - dimethyl acetamide (DMA) is used as the base fluid, (1) and vegetable oil samples at the bottom of the liquid is different, i. e. the operation is not the same conditions,

the existence of more or less error, (2) furthermore, MDEA to use chromatography pure, the price is expensive; otherwise, use the analysis of pure and easy to introduce impurities, increase the complexity of the analysis, (3) boiling point of MDEA, although high, but always more or less volatile, into the chromatographic column, due to the lower column temperature, MDEA residual in the chromatographic column, column efficiency and separation degree, causing unnecessary interference.

One

2 GB in a packed column with DEGS/102. The practice proves that this post is not too high column efficiency, peak shape easy trailing, short service life; therefore Xiangying technology R & D out special packed column xy-6 solvent, solvent peak shape and symmetry, the column efficiency high.

A new method for the improvement of national standard:

(1) without using MDEA, (I) buy city selling pressure oil (solvent-free) as bottom liquid (II) by their own oil as base solution (heated volatile production, no solvent peaks can be, this method, worthy of promotion)
(2) the Tengzhou City Xiangying Analysis Technology Co., Ltd. R & D special packed column xy-6 solvent, solvent peak shape and symmetry, the column efficiency high.

Appendix:

Solvent XY-6 packed column test plant oil example reference spectrum

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Vegetable oil solvent testing for packed column 6

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Before delivery has been aging, use for a period of use after a period of time, such as baseline instability, please aging again.

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Chromatographic column entry termination sampler, the other end without a detector.

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High purity nitrogen, 40 ml/min, flow velocity at the beginning of 50 ° C temperature, with 5 ° C/min to 150 ° C, temperature programmed temperature 12 hours.

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After cooling, the other end received the detector can be used on.

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After aging, should reach the baseline is stable, good separation.

The Chinese technology, analysis technology leader

The gas chromatography technology, the industry benchmark enterprise

1 focused on the development of gas chromatograph

2 committed to becoming a leading supplier of analytical instruments in China

3 to provide you with the most powerful support

4 super high cost quality products

5 high quality and perfect after sales service

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Tengzhou Xiangying Analysis Technology Co. Ltd.

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The phone：0632-5817578, 13963221227

<http://www.xiangying17.net/>