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Determination of formaldehyde and benzene compounds in water by headspace capillary gas chromatography _

In order to allow the masses to drink clean water, the new standard of drinking water has been revised and announced the implementation of the new increase in the 69 toxicological indicators, including benzene series. Benzene can cause cancer, and other benzene compounds are toxic to humans and aquatic plants. Benzene, benzene, toluene, xylene is the most common and most important solvents and chemical raw materials, benzene series in case of pollution of water, it will cause harm to people. Formaldehyde is widely present in the environment, there is a strong stimulus to the human body.

Tengzhou Xiangying Determination Technology Co. Ltd. established benzenes.

Methods: the determination of Benzene Series in water by XiangyingGC7990plus automatic headspace capillary gas chromatography was used to optimize the chromatographic conditions and the headspace conditions.

Results: 8 kinds of components in 7 min to complete the determination and has a good separation effect; in 1 g/L ~ 100 g/L concentration range can get better linearity, correlation coefficient 0.999. The synthetic water samples and standards were determined by the established method, and the variation coefficient was 1.6% ~ 6.5%. The minimum detection concentration is 0.5 g / L (benzene and toluene), 1 g / L (ethylbenzene, styrene, o-xylene and Cumene), 2 g / L (of m-xylene,).

Conclusion: the established method is simple, rapid and has good precision and accuracy. It can be used for the determination of benzene in water

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
Big screen in English two kinds of display, the picture switch is simple and clear, the appearance is fashionable and beautiful.

To improve the automation, intelligent, multi-functional, multidimensional chromatographic system (ARM9-32 bit chip and foreign original software) width of the heating rate, rapid cooling system, temperature control technique for high stability, very good performance price ratio.

Perfect self diagnosis function, can make the user easy to check the fault location and fault type.

Perfect temperature protection and platinum wire resistance to open, short circuit alarm function, to ensure that the temperature is not out of control.

Can be equipped with a built-in AD conversion circuit, can be directly digital output signal to achieve the complete control and analysis of the work on the PC.

Post box can realize negative temperature operation by dry ice or liquid nitrogen.

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