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Determination of the content of B oxygen by Headspace Gas Chromatography

Ethyl cellulose (ethylcellulose) as a pharmaceutical excipient, can be used as drug matrix, coating material, carrier, adhesive and capsule wall material, and so on, the quality standard were received contained in the 2005 edition of the Pharmacopoeia of the people's Republic of China, the United States Pharmacopoeia version 30, 2007 edition of the British Pharmacopoeia and European Pharmacopoeia version 5 [2~5]. The 2005 edition of "Chinese Pharmacopoeia" of the people's Republic of China Annex VII G ethyl cellulose.

Ethyl cellulose (EC) is a thermoplastic nonionic cellulose ether, soluble in various organic solvents, at low temperatures can still retaining the with adhesive, filling and film good flexibility. EC can be used for synthetic plastics, paint, ink, insulating materials such as. Due to lack of toxicity in EC, is applied more and more widely in the medicinal materials and food packaging. Ethoxy in ethylcellulose product content of product.


Tengzhou Xiangying Analysis Technology Co. Ltd. according to the shortages of existing methods for determination of ethoxy ethyl cellulose, a gas chromatography with XiangyingGC7990Plus headspace (HS-GC) determination of ethyl cellulose ethoxyl content method in the vial through hydroiodic acid and ethyl cellulose reaction, the ethoxy ethyl iodide into the reaction. After adding lye and excess hydrogen iodate to prevent corrosion of the injector and the column, followed by headspace gas chromatographic analysis using toluene as internal standard, determination of iodine ethane content by internal standard method, and calculated the ethoxyl content of ethyl cellulose was studied. The hydrolysis temperature and time, headspace sampler in the balance of temperature and time, matrix effect and other factors impact on the determination results showed that the hydrolysis temperature at 140 to 150 DEG C, the temperature changed little, 1.0h can complete hydrolysis reaction ; headspace vial sample in 80 DEG C after 30 min can achieve gas-liquid equilibrium; adding cellulose can eliminate matrix effects and iodine ethane concentration and peak area showed a good linear relationship. In this paper also studied the the method reproducibility and recovery rate, found that the reproducibility of the method better, ethoxy content determination results of the relative standard deviation was 2.83%, recovery rate of 99.6%. This method and the chemical titration method, be consistent with the determination results of samples.

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XiangyingGC7990plus features:

1 excellent performance, competitive price, to meet the customers' choice

2 multi industry adaptability

3 outstanding performance in all kinds of harsh environment

INJ

A variety of sampling ports can be equipped with: packed column inlet, shunt / non shunt capillary inlet,

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