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
Two - butyl ether and chiral impurities by capillary gas chromatography

Dihydroxy butyl ether (DHDB) belongs to duplex promoters of bile secretion, cholagogue, anti-inflammatory, antispasmodic, hepatoprotective, lipid-lowering, a row of stone, increase the secretion of bile effects, and can be used in many kinds of chronic cholelithiasis and cholecystitis, cholangitis, biliary tract infection treatment is a kind of efficient cholagogue, gallstone dissolution agent and liver protecting drugs, it has attracted more and more attention. However, the DHDB structure of the two chiral carbon and potential.


Tengzhou City Xiangying Analysis Technology Co., Ltd. the (CGC) XiangyingGC7990plus capillary column and gas chromatography respectively of dihydroxy butyl ether (DHDB) and chiral impurity separation conditions were investigated and to establish a method for determination of the content of dihydroxy butyl ether (DHDB). (30 m x 0.53 mm x 1.0 m) capillary column chromatography FID detector, internal standard method, DHDB chiral impurity program heating conditions for initial temperature of 120 DEG C, to 0.5 DEG C / min is heated to 180 DEG C, keeping 40 min; DHDB content was measured by temperature programmed conditions to 140 DEG C for 1 min, 20 C / min up to 200 DEG C. DHDB and its isomer can be completely separated; DHDB in 66.70~800.1 $\mu\text{g} / \text{L}$ has a good linear relationship ($y=0.0070x+0.0554$; $R^2 = 0.9996$). The average recovery rate was 99.54% and repeatability RSD was 0.68%, intermediate precision (RSD) was 0.95%, the average content of the samples is 77.92% (RSD=0.25%). The method is simple and rapid determination of chiral impurity analysis and content for DHDB.

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XiangyingGC7990plus features:

- 1 excellent performance, competitive price, to meet the customers' choice
- 2 multi industry adaptability
- 3 outstanding performance in all kinds of harsh environment

INJ

A variety of sampling ports can be equipped with: packed column inlet, shunt / non shunt capillary inlet,

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