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Ammonia, also known as Ammonia water, the main component of $\text{NH}_3 \cdot \text{H}_2\text{O}$, ammonia is an aqueous solution, colorless and transparent and pungent odor. Ammonia -77°C melting point, the boiling point of 36°C , density $0.91 \text{ g} / \text{cm}^3$. Ammonia easily soluble in water, ethanol. Volatile, with a part of the base of the common, ammonia from ammonia into the water obtained. Ammonia toxic, on the eyes, nose, skin irritation and corrosive, can cause suffocation, the air *** high allowable concentration of $30 \text{ mg} / \text{m}^3$. Mainly used as fertilizer.

Industrial ammonia is ammonia 25% to 28% of the aqueous solution, ammonia, only a small part of ammonia molecules react with water to form ammonium ions and hydroxide ions, namely ammonium monohydrate, is only present in the ammonia in the weak base. Ammonia freezing point and concentration of ammonia, the commonly used (wt) 20% concentration of freezing point of about -35°C . Neutralization reaction with acid produces heat. Risk of combustion explosion. The specific heat capacity was $4.3 \times 10^3 \text{ J} / \text{kg} \cdot ^\circ \text{C}$ (10% ammonia).

The determination of industrial ammonia using the traditional chemical acid-base titration, this method is cumbersome, long analysis time, the amount of sample, and the determination of human factors, such as the absorption of gas flow rate, the level of deviation and standard deviation Solution titration, the end of the judgment will result in a large error in the measurement results, and the absorption tube sometimes lax phenomenon of inverted suction, resulting in failure or inaccurate analysis. During the analysis, a large amount of ammonia volatilized is harmful to the health and pollution of the environment.

Some people use ordinary gas chromatographic determination of industrial ammonia, but the ammonia is corrosive, equipment and column life is short, and the peak shape of the tail, quantitative error, no use value, can not promote the use of.

Tengzhou Xiangying Analytical Technology Co., Ltd. is advancing with the times and guided by scientific development. Under the leadership of General Manager Wang Xiaoying, we rely on Xiangying Chromatographic Engineering Research Institute and Analysis and Testing Center, and make unremitting efforts. Out of ammonia for industrial analysis of anti-corrosion gas chromatograph, in a leading domestic level. This instrument to solve the two key core issues: 1. Instrument using anti-corrosion material, especially the thermal conductivity detector

TCD, the use of anti-corrosion silk, enough to resist corrosion of ammonia, life greatly extended. (This technology is the nature of our patented technology, has a very high value); 2. Developed a special anti-corrosion ammonia column, the column-shaped symmetry, column efficiency, long life, is also second to none in the country.

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
www.xiangying17.net

XiangyingGC7990plus features:

1 excellent performance, competitive price, to meet the customers' choice

2 multi industry adaptability

3 outstanding performance in all kinds of harsh environment

INJ

A variety of sampling ports can be equipped with: packed column inlet,
shunt / non shunt capillary inlet,

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